

CardiTrack - MVP 1 User Stories

Extract -- do not edit directly. This file is extracted from the canonical userstories.md; make changes there and re-extract. Release sequencing is governed by the release matrix.

Build status (August 9, 2026): 12 of 17 Figma M1 screens are built (M1-01 through M1-09, plus M1-13 CardiMember Detail, M1-14 Edit CardiMember, M1-15 Device Management); the unbuilt screens are the alert set (M1-10 stub, M1-11/M1-12/M1-16) and M1-17 Health Data Export. Four shipped screens have no Figma M1 frame -- needs design sync: SignInPage, ForgotPasswordPage, VerifyEmailPage, Onboarding/AccountSetupPage (Stories 1.5-1.8).

Stories mapped to MVP 1 screens (M1-01 through M1-17), plus the four shipped screens without Figma M1 frames.

Platform Requirements

- Minimum iOS: 16.0 -- covers ~90%+ of active iPhones; required for modern platform APIs and reliable background push delivery
- Minimum Android: 10 (API 29) -- covers ~85-90% of active Android devices; required for scoped storage and modern permission model
- Target iOS: 18 (latest stable)
- Target Android: 15 / API 35 (latest stable)

Onboarding & Setup

Story 1.1: First-Time User Registration (P0 -- Must Have)

- As a concerned family caregiver
- I want to quickly create an account and understand what CardiTrack does
- So that I can start monitoring my elderly parent's health within minutes
- Acceptance Criteria:
 - Simple signup flow (email/password or social login via Auth0)
 - Clear value proposition on landing page
 - 30-day free trial messaging prominent
 - Account creation leads into email verification (Story 1.7) -- the user cannot enter the app until their email is verified
- Screens: M1-02 (Welcome), M1-03 (Sign Up)

Story 1.2: Adding First CardiMember (P0 -- Must Have)

- As a new CardiTrack user
- I want to easily add my parent as a CardiMember with minimal information
- So that I don't abandon the setup process due to complexity
- Acceptance Criteria:
 - Progressive disclosure (collect basic info first, details later)
 - Required fields: Name, Date of Birth, Relationship
 - Optional fields: Photo, medical notes (encrypted), emergency contacts
 - Clear privacy messaging ("Your parent will be notified")
 - Visual progress indicator (Step 2 of 4)
 - Emergency-phone placeholder localized by device region (PR #8): US/CA "+1 555 000 0000", GB "+44 7700 900000"
-- limitation: all other regions fall back to the US format, notable given the US + EU target market
- Screens: M1-04 (Add First CardiMember)

Story 1.3: Device Connection Wizard (P0 -- Must Have)

- As a caregiver setting up monitoring
- I want to connect my parent's wearable device through a guided wizard
- So that I understand what permissions are needed and why
- Acceptance Criteria:

- Device selection screen with icons (Fitbit -- MVP 1; Apple Watch, Garmin, Samsung -- Coming Soon)
- OAuth flow with clear permission explanations
- "Why we need this" tooltips for each permission
- Success confirmation with sample data preview
- Troubleshooting tips if connection fails
- Support for multiple devices per CardiMember
- Screens: M1-05 (Device Selection), M1-06 (OAuth), M1-07 (Success), M1-08 (Baseline)

Story 1.4: CardiMember Profile Management (P0 -- Must Have)

- As a caregiver
- I want to view and edit a CardiMember's profile (photo, medical notes, emergency contact, monitoring settings)
- So that their information stays accurate and I can quickly act in an emergency
- Acceptance Criteria:
 - View profile summary: name, DOB, relationship, photo, emergency contact
 - Encrypted medical notes (biometric gating deferred to R4 per the release matrix -- not an MVP 1 criterion)
 - Enable/disable monitoring toggle with confirmation
 - Alert sensitivity control (Low / Medium / High)
 - Quick-action buttons: View Dashboard, View Alerts, Manage Devices
 - Danger-zone actions: Pause Monitoring, Remove CardiMember (with confirmation dialogs)
- Screens: M1-13 (CardiMember Detail), M1-14 (Edit CardiMember)

Story 1.5: Sign In (P0 -- shipped; no Figma M1 frame, needs design sync)

- As a returning caregiver
- I want to sign in quickly with my email and password
- So that I can get back to monitoring without friction
- Acceptance Criteria:
 - Email + password form with password show/hide toggle
 - "Remember me" checkbox
 - "Forgot password" link ? Story 1.6
 - Inline error message for failed sign-in (no banner)
 - Social sign-in options (Google / Apple) presented consistently with sign-up
 - Link back to sign-up for users without an account
- Screens: SignInPage (no Figma M1 frame)

Story 1.6: Forgot Password (P0 -- shipped; no Figma M1 frame, needs design sync)

- As a caregiver who forgot my password
- I want to request a reset link by email
- So that I can regain access without contacting support
- Acceptance Criteria:
 - Request state: email input + send-reset-link CTA
 - Confirmation state: "Check your email" with resend option (cooldown to prevent spamming)
 - "Back to sign in" path from both states
- Screens: ForgotPasswordPage (no Figma M1 frame)

Story 1.7: Verify Email (P0 -- shipped; no Figma M1 frame, needs design sync)

- As a newly registered caregiver
- I want to verify my email address right after signing up
- So that my account is secured and I can proceed into the app
- Acceptance Criteria:
 - Hard tenant gate: account creation does not auto-login; the user lands on the verification screen and cannot proceed until verified (threads between Story 1.1 and Story 1.2)
 - "I've verified -- continue" action with a checking state
 - "Open mail app" shortcut
 - Resend verification email (cooldown; confirmation message)

- Clear error state when the address is still unverified
- Screens: VerifyEmailPage (no Figma M1 frame)

Story 1.8: Account-Type Setup (P0 -- shipped; no Figma M1 frame, needs design sync)

- As a first-time user after verification
 - I want to say whether I'm caring for my family or providing care professionally
 - So that CardiTrack can tailor my account
 - Acceptance Criteria:
 - Radio-cards: "My Family" (personal) vs "My Organization" (professional care)
 - Selecting "My Organization" reveals a required Organization Name field
 - Continue disabled until a type is chosen
 - Flagged scope question (not a resolution): the Organization option surfaces business onboarding in MVP 1 while the Guardian Plus business tier is post-R4 in the release matrix -- needs a product decision
 - Screens: Onboarding/AccountSetupPage (no Figma M1 frame)
-

Dashboard & Monitoring

Story 2.1: Daily Health Overview (P0 -- Must Have)

- As a busy family caregiver checking in daily
 - I want to see a quick visual summary of my parent's health status
 - So that I know if everything is okay without reading detailed reports
 - Acceptance Criteria:
 - Traffic light status indicators (Green/Yellow/Orange/Red)
 - Key metrics at-a-glance: Steps, Heart Rate, Sleep Quality
 - "Last synced" timestamp
 - Comparison to baseline ("20% below normal activity")
 - Quick action buttons ("Call Mom", "View Details", "Acknowledge Alert")
 - Screens: M1-09 (Main Dashboard)
-

Alert Management

Story 3.1: Receiving Critical Alerts (P0 -- Must Have)

- As a caregiver receiving an urgent alert
- I want to immediately understand what's wrong and what action to take
- So that I can respond appropriately without panic
- Acceptance Criteria:
 - Alert severity clearly visible (color-coded, icon)
 - Plain language description ("Dad hasn't moved this morning. Typical wake time: 7am. Current time: 11am")
 - Recommended actions ("Call to check in", "Contact emergency services")
 - One-tap actions (Call, SMS, Acknowledge)
 - Alert history visible ("This is the first time this month")
- Screens: M1-10 (Alerts List), M1-11 (Activity Detail), M1-12 (Critical Detail), M1-16 (Heart Rate Detail)

Story 3.3: Alert Acknowledgment & Notes (P1 -- Should Have)

- As a caregiver following up on an alert
- I want to mark it as acknowledged and add notes about my action
- So that other family members know it's been handled
- Acceptance Criteria:
 - Quick acknowledgment button with timestamp
 - Notes field ("Called, he had a cold but is fine")

- Photos upload option (if doctor visit occurred)
 - Alert status: New ? Acknowledged ? Resolved
 - Notification to other family members when acknowledged
 - Screens: M1-10 (Alerts List), M1-11 (Activity Detail), M1-12 (Critical Detail)
-

Settings & Preferences

Story 6.2: Device Management (P1 -- Should Have)

- As a caregiver whose parent switched devices
- I want to disconnect old device and connect new one easily
- So that data continues flowing without interruption
- Acceptance Criteria:
 - List of connected devices with status (Active, Disconnected, Token Expired)
 - Refresh/reconnect button for expired OAuth tokens
 - Primary device designation (when multiple devices connected)
 - Device removal with confirmation ("This will delete connection but keep historical data")
 - Data source indicator on charts (which device provided this data)
- Screens: M1-15 (Device Management)

Story 6.3: Health Data Export (P0 -- Must Have; status: M1-17 not built yet -- entry points show "Coming soon")

- As a caregiver preparing for a doctor's visit or needing records
 - I want to export a CardiMember's health data in standard medical formats
 - So that I can share it with healthcare providers or keep it for my records
 - Acceptance Criteria:
 - Date range selector for the export window
 - Format options: PDF, CSV, FHIR R4 (MVP 1); HL7 v2 (MVP 2); LOINC/CCD (MVP 2); SNOMED CT (MVP 3)
 - Delivery options: save to device, share via system share sheet, email to self
 - Clear format explanations ("FHIR R4 is accepted by most US patient portals and EHR systems")
 - Export confirmation with file size estimate
 - Screens: M1-17 (Health Data Export)
-

Alert Type Stories

Alert Type 1: Activity Alerts (Yellow Severity)

Story 11.1: Gradual Activity Decline

- Display:
 - Chart showing 2-week trend (declining line)
 - Comparison: "Dad's steps: 2,500/day. Normal: 5,000/day (-50%)"
 - Context: "This could indicate illness, pain, or low mood"
- Actions:
 - "Call to check in" (primary button)
 - "Acknowledge" (secondary)
 - "Adjust baseline" (if this is new normal)
- Screen: M1-11 (Alert Detail - Activity)

Alert Type 2: Heart Rate Alerts (Orange Severity)

Story 11.2: Elevated Resting Heart Rate

- Display:
 - Heart rate chart with baseline range shaded

- "Mom's resting HR: 88 bpm. Normal: 68 bpm (+29%)"
- Context: "Elevated for 3 consecutive days. May indicate infection or stress"
- Actions:
 - "Recommend doctor visit" (primary)
 - "Monitor for 2 more days" (secondary)
 - "View detailed HR data"
- Screen: M1-16 (Alert Detail - Heart Rate)

Alert Type 3: Pattern Break (Red Severity)

Story 11.3: No Morning Activity

- Display:
 - Large red alert banner
 - "Dad hasn't moved today. Typical wake time: 7am. Current: 11am"
 - Last known activity timestamp
 - Actions:
 - "Call now" (one-tap phone call, primary)
 - "I'm checking in person"
 - "He told me he'd sleep in today" (dismiss with note)
 - Screen: M1-12 (Alert Detail - Critical)
-

Onboarding Flow UX

Step 1: Value Proposition (30 seconds)

- Hero image: Happy elderly person with Fitbit, smiling family on phone
- Headline: "Peace of Mind for Your Family. \$8/month."
- 3 key benefits with icons:
 - Works with devices they already own (Fitbit, Apple Watch, etc.)
 - AI alerts you BEFORE health issues become emergencies
 - 70% cheaper than medical alert systems (\$8 vs \$47/month)
- CTA: "Start Free 30-Day Trial"

Step 2: Account Creation (1 minute)

- Email/password or "Continue with Google/Apple" (social buttons shipped but not yet wired)
- Checkbox: "I agree to Terms & Privacy Policy" (with links)
- No auto-login after creation -- the user must verify their email (VerifyEmailPage, Story 1.7) before entering the app; PostLoginRouter then routes to account-type setup (Story 1.8) or Add CardiMember

Step 3: Add CardiMember (2 minutes)

- "Who would you like to monitor?"
- Form: Name, DOB, Relationship, Photo (optional)
- Tone: "We'll help you set up monitoring in 4 simple steps" (the wizard is a 4-step flow: Create Account ? Add CardiMember ? Connect Device ? Baseline)

Step 4: Device Connection (3 minutes)

- "What wearable device does [Name] use?"
- Device icons with brands (Fitbit active; others Coming Soon in MVP 1)
- Click ? OAuth flow ? Success
- "Great! We're syncing [Name]'s data. This may take a few minutes."

Step 5: Baseline Learning (Info screen)

- "CardiTrack is learning [Name]'s normal patterns..."
- Progress indicator: "Day 3 of 30"

- "You'll start receiving alerts after we establish a baseline (30 days)"
- Option: "Use statistical alerts in the meantime" (basic threshold alerts)

Step 6: First Dashboard View

- Celebratory tone: "You're all set! Here's [Name]'s health overview."
- The baseline screen's "Invite Family Member First" link ships in MVP 1 but the invite flow (M3-02) is MVP 2 -- the link is currently a dead end
- Guided tour overlay (5 tooltips) -- planned, not shipped (no guided tour exists in the current app):
 - "This shows overall health status"
 - "View detailed trends here"
 - "Alerts appear in this section"
 - "Invite family members here"
 - "Need help? Check our support docs"

MVP 1 Priority Summary

Story	Title	Priority	Build Status
1.1	First-Time User Registration	P0	Built (M1-02, M1-03)
1.2	Adding First CardiMember	P0	Built (M1-04)
1.3	Device Connection Wizard	P0	Built (M1-05-M1-08)
1.4	CardiMember Profile Management	P0	Not built (M1-13, M1-14)
1.5	Sign In	P0	Built -- no Figma M1 frame
1.6	Forgot Password	P0	Built -- no Figma M1 frame
1.7	Verify Email	P0	Built -- no Figma M1 frame
1.8	Account-Type Setup	P0	Built -- no Figma M1 frame
2.1	Daily Health Overview	P0	Built (M1-09)
3.1	Receiving Critical Alerts	P0	Not built (M1-10-M1-12, M1-16)
6.3	Health Data Export	P0	Not built (M1-17 -- P0 but unbuilt)
3.3	Alert Acknowledgment & Notes	P1	Not built
6.2	Device Management	P1	Not built (M1-15)
11.1	Gradual Activity Decline	--	Not built (M1-11)
11.2	Elevated Resting Heart Rate	--	Not built (M1-16)
11.3	No Morning Activity	--	Not built (M1-12)

Telemetry-consent gap (product follow-up): crash reporting and RUM (PR #4) ship enabled with TrackingConsent.Granted hardcoded -- consent is granted by default, with no in-app opt-out and no diagnostics screen. There is no in-app telemetry control in MVP 1. This is in tension with the consent-first principle (Principle 4) -- see Story 7.4 in the canonical userstories.md.

UI/UX Design Principles (MVP 1)

Principle 1: Trust Through Transparency

- Show data source and reasoning for every alert
- Use warm, caring language ("Your mom's activity is lower than usual. Might be worth a check-in call")
- Avoid medical jargon and alarmist language

Principle 2: Simplicity Over Features

- Information hierarchy: Status ? Alert ? Action
- Progressive disclosure (advanced features hidden until needed)

- Mobile-first design

Principle 3: Peace of Mind, Not Panic

- Green status should be prominent when all is well
- Alerts provide context, not just warnings
- Success messaging ("Your dad had his most active week this month!")

Principle 4: Respect for Elderly Dignity

- Never use patronizing language or imagery
- Focus on independence and wellness, not decline
- Consent-first approach to all monitoring

Principle 5: Multi-Generational Accessibility

- WCAG AA compliance minimum (AAA preferred)
- Font size options (small/medium/large)
- High contrast mode
- Keyboard navigation support
- Screen reader optimization

Source: Extracted from userstories.md v1.1 (manually re-synced August 9, 2026)

Screens covered: M1-01 through M1-17, plus SignInPage, ForgotPasswordPage, VerifyEmailPage, and Onboarding/AccountSetupPage (no Figma M1 frames -- need design sync)